

Digital Analogs for Certain[®] and External Hex Connections



World Class **Flexibility**

The New **Digital Analog** Offers A Unique One Stop Solution For Your Various Digital Workflow Needs

1 Unparalleled Precision

The Certain and External Hex Digital Analogs are designed and manufactured to exact tolerances and specifications in order to accurately replicate implant placement in a 3D printed model and stone model.

2 Dual Functioning Design

Our re-designed Digital Analogs can now be used in both a stone model or a 3D printed model workflow for cases using Certain and External Hex Implants.


Note: Now Available in a Carbon printed 3D model

 **ZIMMER BIOMET**
Your progress. Our promise.®

Digital Analogs for Certain and External Hex Connections


NEW DESIGN

The new design incorporates the Zimmer Biomet Logo for easy identification of an OEM (Original Equipment Manufacture) product.



CERTAIN CONNECTION

The Certain Analog internal connection is color-coded by size for easy identification in the lab and the model.



EXTERNAL HEX CONNECTION

Additional step in the apical design of the analog indicates the newer design that can be used in a 3D printed model.

SAME PART NUMBERS

No changes to part numbers, same products you have been ordering for years but now with greater functionality.

Certain Internal Connection Digital Analogs

Seating Surface	Item No.
3.4 MM(D)	IMMILA
4.1 MM(D)	IILA20
5.0 MM(D)	IILAW5
6.0 MM(D)	IILAW6

External Hex Connection Digital Analogs

Seating Surface	Item No.
3.4 MM(D)	MMILA
4.1 MM (D)	ILA20
5.0 MM(D)	ILAW5
6.0 MM(D)	ILAW6



Contact us at 1-800-342-5454 or visit zimmerbiometdental.com

Zimmer Biomet Dental
Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Phone: +1-561-776-6700
Fax: +1-561-776-1272

Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured by one or more of the dental subsidiaries of Zimmer Biomet Holdings, Inc. and marketed and distributed by Zimmer Biomet Dental and its authorized marketing partners. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental. ZB0709 REV B 09/19 ©2019 Zimmer Biomet. All rights reserved.

